



The Mediwet humidifiers are developed for oxygen therapy at hospital or at home. This type of equipment is indispensable to prevent dehydration of the mucous membranes when inhaling oxygen during an extended period of time, or at large flows. It is intended to be connected to the outlet of a wall flowmeter, a regulator flowmeter, or an oxygen concentrator.

Humidifiers Mediwet

EDITION 2/2012

The reusable humidifier is ergonomic, comfortable, economic, and easy to maintain:

- Reduced dimensions and weight.
- Silent as a result to a high performance diffuser.
- All components, except the diffuser, are reusable for minimum 20 autoclave cycles.
- Improved adaptation to the needs and possibilities of the user with two autoclave possibilities, depending on the material.
- Fast to assemble/disassemble.

MEDIWET 2 HUMIDIFIER

Art. No.	Denomination
149106	Mediwet 2 134°C inlet connection G3/8
149107	Mediwet 2 134°C inlet connection UNF 9/16 F
150552	Bottle gasket 134 °C, 10 pcs



Technical Data - MEDIWET 2 HUMIDIFIER

Contents:	Only sterile water
Dimensions:	
Height:	190×67 Width (incl. hose nipple) ×57Ø mm
Weight:	Polysulphone version 115 g
Capacity:	200 ml of water
Consumption:	8 ml of water per hour at a as flow of 10 l/min at 20 °C
Material:	Jar autoclavable at 134 °C: polysulphone Lid and outlet hose nipple: polypropylene Inlet nut: Chromed brass Diffuser: Polyethylene Hose, flat gasket: Silicone Toric gasket: Ethylene propylene
Outlet connection	Tapered hose nipple for hose 6×9mm (recommended length 2m)
Cleaning:	Water, non abrasive detergent. Never use solvents.
Disinfection:	An alcoholic solution, or other solution compatible with the material according to the disinfectant manufacturer.
Autoclave:	Polysulphone version 134 °C for 18 min.
Special case:	Diffuser: exchange at every cleaning - do not sterilise!
Maintenance	Check that the humidifier is whole and air tight before use. Monthly exchange of gaskets is recommended. Exchange the diffuser when its microperforations no longer exist.
Durability	Minimum 20 autoclave cycles under the condition that all instructions accompanying the product are adhered.

Recycling

The product shall be recycled according to local regulations.



Air Liquide Gas AB
Box 2911
212 09 Malmö
Sverige
Tel: +46 40 38 10 00
Fax: +46 40 43 69 43
www.airliquide.se

Air Liquide Danmark A/S
Høje Taastrupvej 42
2630 Taastrup
Danmark
Tel: +45 43 55 50 50
Fax: +45 43 55 50 65
www.airliquide.dk

Air Liquide Norway AS
Ryghgata 2 B; Postboks 243
3051 Mjøndalen
Norge
Tel: +47 32 27 41 40
Fax: +47 32 27 41 59
www.airliquide.no